

## Molecular cloning methods in our hands (Version 1.0)

### 1. Materials and reagents

High-Fidelity PCR polymerase: 2× Hieff Canace® PCR Master Mix (Yeasten, Cat# **10136ES01**)

High-Fidelity PCR polymerase: PrimeSTAR max 2× Master Mix (Takara, Cat# **R045Q**)

Colony PCR polymerase: 2x Es Taq MasterMix (CWBIO, Cat# **CW0690**)

DNA ligase: T4 DNA Ligase (NEB, Cat# **M0202M**)

DNA purification (Small-scale): BBI DNA purification Kit

DNA purification (Large-scale): 96-Well Plate PCR Purification Kit (BBI, Cat# **B518145**)

DNA gel recovery: QIAquick Gel Extraction Kit (QIAGEN, Cat# **28706**)

Fragments recombination: Hieff Clone® Plus One Step Cloning Kit (Yeasten, Cat# **10911ES20**)

*# You may need to read through the commercial protocol before cloning.*

### 2. Gene isolation via PCR method

#### ➤ PCR reaction Mix

■ Template DNA (1 ng/μL) <sup>#1</sup>	2 μL
■ FP (10 μM)	2 μL
■ RP (10 μM)	2 μL
■ DMSO <sup>#2</sup>	4 μL
■ 2x Canace Mix <sup>#3</sup>	50 μL
■ PCR grade water	to 100 μL <sup>#4</sup>

*#1 Plasmid DNA should be diluted to <1 ng/μL, high concentration of template may inhibit PCR reaction.*

*#2 DMSO is a frequently-used PCR additive. 1-3% of DMSO may increase the specificity or may help you to deal with high GC-contained templates.*

*#3 We have 2×ready-to-use master mix which contains polymerase, dNTPs, Buffers. Or you need to add 1 U/50 μL of Canace polymerase additionally.*

*#4 Commonly-used thermo cyclers are applicable for 50μL per hole, **EXCESS IS FORBIDDEN**.*

#### ➤ PCR cycling condition<sup>#4</sup>

■ 98°C	3min	} 31-35 cycles
■ 98°C	10sec <sup>#2</sup>	
■ 60°C <sup>#1</sup>	20sec	
■ 72°C	30sec/kb <sup>#3</sup>	
■ 72°C	5min	

*#1 The optimization of annealing temperature should be started from 60°C.*

*#2 Longer denaturing time may help you to deal with high GC-contained templates.*

*#3 For simple, low GC, plasmid DNA, the extension time can be decrease as 10 sec/kb.*

*#4 For **PrimeSTAR** max, we have another well-established protocol to use.*

### 3. Plasmid linearization via restriction digestion

#### ➤ Digestion reaction Mix

- 
- |                       |                |
|-----------------------|----------------|
| ■ DNA                 | 1 µg           |
| ■ 10×Buffer           | 50 µL          |
| ■ Enzyme              | 1 µL (or 10 U) |
| ■ Nuclease-free water | To 50 µL       |

- Incubate the reaction mix at suitable temperature.

# *5–10 units of enzyme per µg DNA in a 1 hour digest. Overnight incubation may need to guarantee totally digesting*

#### 4. DNA electrophoresis

- Mix 1 µL of loading dye with 4 µL of DNA sample for electrophoresis
- Gel preparation.
- Prepare gel mould and combs
  - Boiling gel, add **10000×** of gel dye, mix well
  - Pooling gel, cooling for 30 min

# *1% of agarose is suitable for detecting DNA fragments >100bp*

#### 5. DNA purification

- In a 0.2 mL PCR tube, add 5 volumes of **Binding Buffer** (with 25% of isopropanol) to each volume of DNA sample.
- Mix well by pipet, transfer sample mixtures to the wells of 96-well **DNA purifier**.
- Centrifuge at  $\geq 3,000 \times g$  (4000 rpm) for 5 minutes until sample mixtures have been completely filtered. Discard the flow-through.
- Add 750 µL **Wash Buffer** to each well. Centrifuge at  $\geq 3,000 \times g$  for 5 minutes.
- Using air-dryer to dry the columns.
- Add 30-40 µL pre-heated water or **Elution Buffer** directly to the column matrix in each well. Incubate the plate in the 55°C incubator.
- Transfer the DNA purifier onto an Elution Plate and centrifuge at  $\geq 3,000 \times g$  for 3 min to elute the DNA.

# *This protocol is modified for large-scale assay (>10 samples)*

# *For small-scale assays, please refer to kit manual*

#### 6. DNA ligation

- Carefully DNA cleaning-up should be performed to avoid high false positive. Gel purifying the backbone and insert DNA if it is possible.
- Mixture condition (backbone : insert = 1:3, for 2 fragments rec.)
- Backbone:  $[0.02 \times \text{basepairs}]$  ng
  - Insert:  $[0.06 \times \text{basepairs}]$  ng
  - Add up to 4 µL, mix with 4 µL of 2 × Enzyme mix.

# *Normally, we use 40 ng of backbone for recombination.*

#### 7. Transformation

- Take 1 µL of ligation mixture from above, add 10 µL of competent cells into each well.
- Mix well and incubate on ice for 20 min. Thermal shock under 42°C for 1 min, incubate on ice for 2 min. Add 160 µL of LB (without antibiotics), incubate for 1 h under 37°C.
- Take 20 µL<sup>#1</sup> of cell suspensions for inoculation.

#1 *For huge constructs like Huc, GFAP, inoculating more.*